



VETERINARY FLUID WARMER Q-Warm

FEATURES

Open system:

- To use independently or be equipped to other warming system, infusion pumps or feeding pump
- Accepts standard IV tube, no special disposables needed

Heating up to the patient:

• IV tube is completely wrapped in, no heat loss

One-touch-control:

- Temperature automatically controlled
- Easy to operate
- Effective and gentle heating
- Compact in size

TECHNICAL SPECIFICATIONS Model: Q-Warm

Wodel. Q-Walli

Temperature setting: 39°C ex-works

Power supply: a.c.100~240V/50~60Hz

Power consumption: Max. 50VA

Type of protection against electric shock: Class I

Degree of protection against electric shock: BF Applied part; Defibrillation-protected

Degree of protection against ingress of liquids: IPX2

Temperature accuracy: ±2.0°C

Overheat protection: 42°C

Low temperature alarm: 37°C

Warming up time: From 20°C to 39°C approx.3 min.

Operating mode: Continuous

Dimension (W*D*H): 70×25×90mm

Net Weight: 0.7kg

Warming profile: length 1000mm; compatible to IV tube 3.5-5.0mm O.D.

* Subject to technical modifications without notice





