### Transducers C5-2Q MC8-4Q MC9-3TQ C6-2MQ 2 - 5 MHz 4 - 8 MHz 3 - 9 MHz 2 - 6 MHz 3D/4D Mechanical Curved Array Micro - Convex Array Micro - convex Array Curved Array P7-3Q L12-5Q L17-7HQ P5-1Q 7 - 17 MHz 3 - 7 MHz 1 - 5 MHz 5 -12 MHz

Phased Array

# Flexible and Available at Fingertips

Acclarix LX3 VET Diagnostic Ultrasound System



Linear Array

### Global Headquarters:

**High Frequency** 

Linear Array

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126 +1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability. ENG-US-LX3 VET-V1.0-20210818

Pediatric Phased Array







## Stunning Clarity

Compacted with innovative Edan TAI technology and multiple imaging processing technologies, Acclarix LX3 VET could perfectly display ultrasound image in different modes, assisting sonographer to make more precise diagnosis.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
- In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality





Canine Heart (M mode)

### Canine Heart (B mode)





Gall bladder (B mode)

### Distinctive Design





design for better experience

user panel to simplify diagnosis

### New Generation of Capability and Flexibility

Integrated with compact appearance, cutting-edge technology and intelligent workflow, Acclarix LX3 VET is to provide more valuable innovation. Its excellent image performance plays an important role in the prevention and early detection in primary medical care, enriching the diagnostic application and enhancing the clinical value.

### The Virtue of Value

The Acclarix LX3 VET delivers unmatched value and performance across a broad range of applications:

- Canine
- Feline
- Bovine Abdomen
- Ovine
- Equine

- Reproduction
- Cardiac


- Superficial



- In Doppler, TAI automatically adjusts for flow state, providing improved continuity, border detection and filling of color



Liver (C mode)



kidney (B mode)

Kidney (C mode)





5 active tra

