



URIT-880

Chemistry Analyzer

- High performance
- Easy touch screen operation
- Convenient maintenance design
- Support barcode reader
- With 8 inner incubation positions

URIT

URIT-880


Chemistry Analyzer

| Specifications |

Analysis method:	Kinetic, end point, two point, dual wavelength, ABS
Light source:	halogen lamp, 6V/10W
Wavelength (in nm):	340, 405, 492, 510, 546, 578, 630, 700
Absorbance range:	-0.3~4.0Abs
Wavelength accuracy:	± 2nm
Spectrum band:	± 10nm
Flowcell:	10mm quartz cell, 32 μ L volume
Temperature control:	room temperature, 37°C ($\pm 0.1^\circ\text{C}$ variable)
Sampling volume:	100 μ L ~ 9999 μ L adjustable, recommended: 500 μ L
Cross contamination:	$\leq 1.0\%$
Stability:	$\leq 0.005\text{Abs/h}$
Memory:	Up to 300 items and 30000 test results
Display:	Large LCD touch screen display
Printer:	Built-in thermal printer
Interface:	Parallel port (For external printer), RJ-45, 3USB, SDID, PS-2
Ambient:	Room temperature: 10 $^\circ\text{C}$ ~ 32 $^\circ\text{C}$; Humidity: $\leq 40\%$ ~ 85%
Power:	AC (110-220)V, 50Hz
Dimension:	362.5 x 341 x 165mm (length * width * height)



7 inch touch screen



8 inner incubation



Multi function ports

| Features |

- High-quality grating & close optical system
- Easy touch screen operation with 7 inch colour screen
- Convenient maintenance design
- Large memory capacity
- With 10 inner incubation positions
- Support barcode reader

URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Head Office

Add.: No.4 East Alley, Jiuhua Road, Guilin,
Guangxi 541001, P.R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

E-mail: export@uritest.com service@uritest.com

Shenzhen Office

Add.: Room E.J, 24th Floor, New Baohui Building, Nanhai Road,
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

Http://www.urit.com

