




VETCORDER™

EVERY PATIENT | EVERY TIME




Every Patient. Every Time.

Improve the quality of patient care with a portable monitoring tool smaller than a credit card. VetCorder™ is a powerful, state-of-the-art multi-parameter patient monitor that conveniently and accurately tracks SpO2, heart rate and ECG to help ensure patient safety from admission to discharge.

A modern monitoring method that's simple, flexible and practical, VetCorder is both convenient and cost effective. It can be used during routine exams, emergencies, advanced and elective surgeries, dental treatments, recovery or any other procedures requiring anesthesia. Using VetCorder is simply the right thing to do for every patient. Every time.





Shown
actual size

By providing enhanced monitoring capabilities, Vetcorder improves patient safety and instills confidence in staff and animal owners. The pocket-sized touchscreen instrument, featuring Bluetooth connectivity to iOS and Android mobile devices, does the job of multiple large monitors often exclusive to surgical suites. Affordable and economical at just \$995 (MSRP), Vetcorder can help both new and established veterinary practices looking to optimize monitoring protocols.

Monitoring Functions:

- SpO2
- Heart rate
- ECG
- Ability to provide convenient surface temperature spot check
- ECG and Pleth Waveforms

Product Features:

- At a 3.5" x 2.2" x .5" size, it fits in the palm of your hand (and easily in your pocket)
- Weighs under 80g – less than your typical smartphone!
- Bluetooth-enabled to allow continuous remote monitoring with your smartphone or tablet up to 30 feet away
- Battery longevity for over 12 hours of continuous monitoring and can hold its charge for up to 3 months on standby



Uses:

Vetcorder can be used on a wide variety of animal species, including dogs, cats, other companion animals, livestock, non-human primates, reptiles, birds and exotics. It is ideal for:



Specifications Chart:

SpO2	
Range	70% - 100%
Accuracy	80-100% +/- 2%, 70-79 +/- 3%
HR/PR	
Range	30 to 350 bpm
Heart Rate	+/- 2 bpm or +/- 2%
Connectivity	
Wireless	Bluetooth 4.0
Battery	
ECG Recorder	Over 1,000 times
Bluetooth On	12+ hours
Battery Type	Rechargeable lithium-polymer 560 mAh

Each System Includes:

- One Vetcorder
- SpO2 lingual sensor
- ECG cable with two sets of clips
- Multi-parameter adapter
- Silicon protective boot
- USB charging cable with AC adapter
- Quick-start guide
- Warranty
- In-service training

Using Vetcorder from admission to discharge helps veterinary practices achieve a new level of patient security and peace-of-mind. The worst message to deliver to a client is the loss of a loved companion, especially when it could have been avoided. Vetcorder keeps staff informed and connected to vital information to help ensure patients are not compromised.

Bluetooth connectivity via mobile device allows veterinarians, technicians and staff the flexibility to monitor patient vitals while executing other tasks or working with other patients in the practice or hospital.

Vetcorder is widely used in:

- Veterinary practices and hospitals
- Exotic pet centers
- Emergency animal facilities
- Laboratories
- Zoos
- Universities
- Production animal environments

World Class Customer Support:

Every Vetcorder system is backed by Sentier's dedicated service and support team. Purchase includes free in-service training with one of our technical experts and access to technical support via phone or web site. The Sentier team is always one call away to assist your practice or facility with questions, uses or technical support.

Bottom Line Patient Security:

Vetcorder pays for itself in a short amount of time via minor monitoring fees. It is a rare occurrence that a pet parent would not want their loved family member both cared for and monitored from exam to recovery. It is the right thing to do for patient care.



The Case for Vetcorder in Recovery:

According to a recent study, more than half of the deaths tied to veterinary surgery happened within three hours of recovering from anesthesia, which is used in a multitude of minor procedures.¹

In a 2016 study, 57 percent of dogs and 35.5 percent of cats that underwent procedures using anesthesia experienced an adverse event, including: arousal or breakthrough pain (14.9%), hypoventilation (13.5%), hypotension (10.3%), arrhythmias (5.8%), hyperthermia/hypothermia (5.0%), airway complications (4.8%), recovery excitation (4.6%), aspiration risk (4.5%), desaturation (2.8%), hypertension (1.7%) and other ailments (3.7%).²

The study found that using monitoring tools like a Vetcorder can be easily integrated in vet practices to provide valuable evidence and data for veterinarians, veterinary technicians and anesthesiologists to maintain optimal health for their patients during procedures and in recovery.

¹Brodbeck, DC et al (2008) "The Risk of Death: The Confidential Enquiry into Postoperative Small Animal Fatalities." *Veterinary Anesthesia and Analgesia* 35, 365-373.

²McMillan, M et al (2016) "Adverse event surveillance in small animal anesthesia: an intervention-based, voluntary reporting audit." *Veterinary Anesthesia and Analgesia* 43(2): 128-35.



Why Vetcorder?

Industry professionals say it best.

“Vetcorder does something no other monitor can do for us. It shows that status of our patients in recovery. Every other solution we’ve tried has fallen short of Vetcorder’s versatility and value. **With Vetcorder, our clinic can deliver a higher quality of patient care.**”

-Taylor Stewart, RVT

“After a week of using it, my team of **8 CVTs asked me to buy more units since they were fighting over the one I bought.** We now have a total of 3 and my other machines are gathering dust. The Sentier team has also been great about answering my questions that we have. My techs and I highly recommend the Vetcorder.”

-Dr. Ilan Waskow, Mile High Veterinary Hospital

“Continuous monitoring is especially important when treating with drugs that affect the heart. Thanks to the Vetcorder, **I have peace of mind during transport and imaging.**”

-Shalee Ratliff, RVT

“We just can’t get enough of this device! **Use it in the ER, out in the field, post operatively or even during the intraoperative period this device is incredibly versatile.** Did we mention its Bluetooth feature?”

-Veterinary Anesthesia Nerds “Top 10 Favorite Things” 2016, 2017 and 2018!



Do the right thing for patients and clients.

Leverage the flexible, convenient power of Vetcorder for admission to discharge patient monitoring, ensuring animal safety and well-being. With Vetcorder, it’s easy to monitor.

Every patient. Every time!

To learn more, call your preferred distributor and set up a free two-week trial today or contact us at:
1-844-VETCORDER (1-844-838-2673)
info@sentierconnect.com
www.sentierconnect.com

©2019 MAI Animal Health. All rights reserved.



VETCORDER™
EVERY PATIENT | EVERY TIME

 **SENTIER**

 **MAI** ANIMAL HEALTH